

Quantum Imaging and Therapeutic Associates

Ankle – 1.5 Tesla

Indications	OCD, tendon pathology, occult fracture, mass
Coil	Ankle, flex L
Position	Supine, ankle at 90 degrees
FOV	14-16 cm
Comments / Special	Skin marker at site of pain, palpable lesion, scar

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax PD	3000	20-40			3/1	512 x 256
Ax T2 FS	>2000	50-70			3/1	256 x 256
Cor PD FS	3000	20-40			3/1	256 x 256
Sag T1	400-800	Min.			3/1	256 x 192
Sag STIR	>2000	20-40	150		3/1	256 x 192