

Quantum Imaging and Therapeutic Associates

Brain MRA – 1.5 Tesla

Indications	Cerebral aneurysm, vascular stenosis, thromboembolic disease.
Coil	Head
Position	Supine
FOV	22 cm
Comments / Special	MIP 12 projections @ 15 degree increments right to left and anterior to posterior. MIP left anterior, right anterior, and posterior circulations 12 projections @ 15 degree increments right to left.

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax 3D TOF						