

Quantum Imaging and Therapeutic Associates

Elbow (Arthrogram) – 1.5 Tesla

Indications	Ligament evaluation, intra-articular loose body
Coil	Knee, flex M
Position	Supine, arm extended above head, palm down
FOV	12-14 cm
Comments / Special	N/A

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax T1	400-800	Min.			3/1	256 x 256
Ax T2 FS	>2000	50-60			3/1	256 x 256
Cor T1	400-800	Min.			3/1	256 x 256
Cor T2 FS	>2000	50-60			3/1	256 x 256
Sag T1 FS	400-800	Min.			3/1	256 x 256
Sag 3D GRE	30-40	Min. (in phase)		30	1.5/0	256 x 256