

# Quantum Imaging and Therapeutic Associates

## Female Pelvis – 1.5 Tesla

Indications	Pelvic pain, uterine or adnexal disease, uterine anomaly.
Coil	Flex L, torso
Position	Supine
FOV	32 cm (All T1, Cor T2), 24 cm (Sag T2, Ax T2)
Comments / Special	Cor T2 should include bottom of kidneys. If specifically for uterine anomaly, known or suspected, axes should be oriented to the uterus.

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Sag T2	>6000	100-120			5/1	320 x 320
Ax T2	>6000	100-120			5/1	320 x 320
Cor T2	>6000	100-120			5/1	256 x 224
Ax T1	500-700	10-20			5/1	256 x 256
Ax T1 FS SPGR	200	1.5			5/1	256 x 256
Ax T1 FS SPGR gad	200	1.5			5/1	256 x 256
Sag T1 FS SPGR gad	200	1.5			5/1	256 x 224
Cor T1 FS SPGR gad	200	1.5			5/1	256 x 224