

# Quantum Imaging and Therapeutic Associates

## Liver Transplant – 1.5 Tesla

Indications	Liver transplant evaluation.
Coil	Flex, torso
Position	Supine
FOV	32 cm
Comments / Special	Do not use Eovist. Must use contrast agent with non-biliary excretion.

<u>Sequence</u>	<u>TR</u>	<u>TE</u>	<u>TI</u>	<u>Flip</u>	<u>ST/sp</u>	<u>Matrix</u>
Ax dual echo T1 in/out phase						
Cor T2 SSFSE						
Ax T2 SSFSE						
Ax T2 SSFSE FS						
Ax T1 FS SPGR						
Ax T1 FS SPGR post gad (20s, 60s, 120s)						
Cor T1 FS SPGR post gad (180s)						